


Switch Accessible Apps for iPad/iPhone




Jane Farrall, Jane Farrall Consulting www.janefarrall.com
Alex Dunn, Inclusioneers <http://smartinclusion.wikispaces.com/>
9th June 2013




We have divided this list into two sections. The first is a listing of switch accessible Apps that aren't primarily designed for AAC. The second section is a listing of switch accessible Apps which are primarily designed for AAC. The only App that is in both sections is Scene&Heard. Although Scene&Heard is designed primarily as an AAC App, we are both using it for a range of educational reasons and have therefore added it to both lists.




This list is updated regularly and located at <http://www.janefarrall.com/html/resources/Switch%20Accessible%20Apps%20for%20iPad.pdf>




Switch Accessible Apps (non-AAC)




App	Device	About the App	Keystrokes	Switch Access
 A Leader Is Free	iPhone, iPad	A talking book about leaders and leadership. SymbolStix symbols illustrate key words. Users can also play bingo with symbols in the text, match symbols and play concentration with them. Part of the Click to Read series.	Switches 1 – 4 can be any number key, ~ + 1, space or enter.	1 switch scanning. 2 or 3 switch step scanning. Auditory scanning option. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.




App	Device	About the App	Keystrokes	Switch Access
 <p>Adapted Play Book – Dirty Duds \$5.99 (CA/US) \$6.49 (AU) £3.99 (UK)</p>	iPad	An interactive play book based around the Dirty Duds book by Pati King DeBaun. Four activities are included. These are: a simple book which introduces us to the story; an activity where the user has to find the correct picture to match a word from the story from a choice of 6; an errorless activity to explore the story items; a song activity based on the book where the user can pick one of 6 items to have a verse sung.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 <p>Adapted Play Book – EIEIO! \$5.99 (CA/US) \$6.49 (AU) £3.99 (UK)</p>	iPad	An interactive play book based around “Old McDonald Has a Farm”. Four activities are included. These are: a simple book which introduces us to the animal sounds; an animal sounds activity where the user has to find the correct sound for each animal from a choice of 6; an errorless activity to explore the animal sounds; a song activity with the song “On the farm, look what I see”, where the user can pick one of 6 items to have a verse sung.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 <p>Adapted Play Book – Hurry! Hurry! \$5.99 (CA/US) \$6.49 (AU) £3.99 (UK)</p>	iPad	An interactive play book around peanut butter. Four activities are included. These are: an engaging and simple book which talks about the character (a boy) and where he likes peanut butter e.g. on his face; a language activity where the user has to find the picture that goes with a spoken word from a choice of 6; an errorless activity to explore the places the character likes peanut butter ; a song activity with the song “I like peanut butter”, where the user can pick one of 6 items to have a verse sung.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.




App	Device	About the App	Keystrokes	Switch Access
 Adapted Play Book – Peanut Butter \$5.99 (CA/US) \$6.49 (AU) £3.99 (UK)	iPad	An interactive play book around peanut butter. Four activities are included. These are: an engaging and simple book which talks about the character (a boy) and where he likes peanut butter e.g. on his face; a language activity where the user has to find the picture that goes with a spoken word from a choice of 6; an errorless activity to explore the places the character likes peanut butter ; a song activity with the song "I like peanut butter", where the user can pick one of 6 items to have a verse sung.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 Aunty Maggie's Recipe \$2.99 (CA/AU/US) £1.99 (UK)	iPad	This app from Inclusive Technology has potion brewing as its goal! Small targets and multiple choices make it ideal for students who need to learn to make choices and to practice touch accuracy or scanning.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	1 and 2 switch visual scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 Big Bang Patterns \$13.99 (CA/US) \$14.99 (AU) £9.99 (UK)	iPad	Simple cause and effect activities, specifically designed to provide for people with complex needs and visual impairments. The graphics are simple and based on geometric shapes. When the user presses the switch, the graphic is animated and is accompanied by music. The background colour and the colour of the graphics can be changed. 4 different sets of graphics are included within the app.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	1 and 2 switch access. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.




App	Device	About the App	Keystrokes	Switch Access
 Big Bang Pictures \$13.99 (CA/US) \$14.99 (AU) £9.99 (UK)	iPad	Simple cause and effect activities, specifically designed to provide for people with complex needs and visual impairments. The graphics are simple and include images such as a bus, sun, person running. When the user presses the switch, the graphic is animated and is accompanied by music. The background colour and the colour of the graphics can be changed. Two extra activities include the same images, but presented as line drawings in either a simple or more complex format. 1 or 2 switch scanning.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	1 and 2 switch access. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 Big Trucks \$2.99 (CA/AU/US) £1.99 (UK)	iPad	A simple one or two switch game. Two “big trucks” e.g. digger, fire truck appear onscreen. Switch 1 causes the truck on the left to animate (and occasionally be replaced). Switch 2 causes the truck on the right to animate (and occasionally be replaced).	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	1 or 2 switch cause and effect. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 Bumper Cars \$2.99 (CA/AU/US) £1.99 (UK)	iPad	A simple game for one or two players – although it works best with two. Each player has a bumper car. During their turn they need to hit their switch or touch the bumper car on the screen. The cars then bump and the player is assigned a point value by the app. First player to 10 wins. User has no control over the performance of the bumper car, hence the game is described as errorless.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Visual scanning with one switch per player. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.



App	Device	About the App	Keystrokes	Switch Access
 <p>Catch the Cow \$1.99 (CA/AU/US) £1.49 (UK)</p>	iPad	Catch the Cow aims to help users learn to select pictures using a scanning device. It aims to provide teachers, clinicians or parents with an App to assess beginning skills for choice making through auto scanning. Direct touch can also be used and accuracy of the two modes can be compared.	Switch – Space or ~ + 1	Offers 1 switch visual scanning. Can be accessed with RJ Cooper/Therapy Box, Blue2 & APPLICATOR switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.
 <p>Cause and Effect Sensory Light Box \$2.99 (CA/AU/US) £1.99 (UK)</p>	iPhone, iPad	A simple cause and effect app with a number of visual effects that occur when the user presses their switch. Offers a choice of 21 different visual patterns. App can also be activated by sound (using the microphone) or by the screen.	Switch – Space or ~ + 1	Offers 1 switch visual access. Can be accessed with RJ Cooper/Therapy Box, Blue2 & APPLICATOR switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 <p>Five Little Aliens \$2.99 (CA/AU/US) £1.99 (UK)</p>	iPad	A simple to play activity app from Inclusive Technology that requires just a touch anywhere on the screen or a switch press to continue the song. Ideal for young children and those requiring simple access or learning numbers 1 - 5. The music section and activities are cause and effect.	Switch – space or ~ + 1	Offers 1 switch visual access. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.




App	Device	About the App	Keystrokes	Switch Access
 <p>Five Sharks Swimming \$2.99 (CA/AU/US) £1.99 (UK)</p>	iPad	A simple to play activity app from Inclusive Technology that requires just a touch anywhere on the screen or a switch press to continue the song. Ideal for young children and those requiring simple access or learning numbers 1 - 5. The music section and early activities are cause and effect, with the final activity introducing some 1 or 2 switch visual scanning.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	1 and 2 switch visual scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 <p>Hurdle Champion \$2.99 (CA/AU/US) £1.99 (UK)</p>	iPad	A simple single switch game with some timing requirements. Can be played at easy, medium or hard levels. Players need to press their switch within a given time frame to get the hurdler to jump the hurdle. In the easy level, the time frame within which they need to press the switch goes for a couple of seconds, but if you increase the level of difficulty then this time frame decreases.	Switch – space or ~ + 1	Offers 1 switch visual access. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 <p>I Can Help \$9.99 (CA/US) \$10.49 (AU) £6.99 (UK)</p>	iPhone, iPad	A talking book about helping. Key words are illustrated by SymbolStix. Users can also play bingo with symbols in the text, match symbols and play concentration with them. Part of the Click to Read series.	Switches 1 – 4 can be any number key, ~ + 1, space or enter.	Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.




App	Device	About the App	Keystrokes	Switch Access
 <p>Jungle Adventure \$2.99 (CA/AU/US) £1.99 (UK)</p>	iPad	An arcade style single switch game with timing requirements. The user needs to press the switch at the right time to complete each activity and collect the 6 golden coins. The activities include jumping a hungry crocodile, leaping a rolling log and swinging over a chasm.	Switch 1 – space or ~ + 1	Offers 1 switch visual access. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 <p>Koppy Kattz \$19.99 (CA/US) \$20.99 (AU) £13.99 (UK) (Lite version available for \$0.99)</p>	iPhone, iPad	Based on the old Simon Says game – four coloured buttons appear on the screen. In cause and effect mode, the user can select anyone of the four buttons to have the colour highlight and the associated tone play. In Monkey See mode the user sees a sequence play and then copies it. The user is given multiple tries at replicating just one selection, before moving onto a sequence of two. They are then give multiple practices at replicating a two button sequence before moving onto a 3 button sequence, etc., through to a 5 button sequence. If they get a sequence wrong, then the sequence is shortened by 1 and they practice again at the earlier level. In Koppy Katt mode, a sequence of 1 to 5 buttons is played and the user must replicate it to stay in the game. Multi-player mode is the same as Koppy Katt mode but with multiple players. In the settings, the length of the starting sequence in Monkey See and Koppy Katz can be controlled. Default difficulty levels of easy, medium and hard can also be changed in the settings.	Switches 1 – 4 can be any number key, ~ + 1, space or enter.	1, 2, or 3 switch visual scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 <p>Little Lost Penguin \$2.99 (CA/AU/US) £1.99 (UK)</p>	iPad	An arcade style single switch game with timing requirements. The user needs to press the switch at the right time for the penguin to get past obstacles such as an angry polar bear and a moving ice shelf. Along the way, the penguin collects fish to sustain him. At the end, the penguin finds his way home to Mum and Dad and the penguin clan who do a welcome home dance.	Switch 1 – space or ~ + 1	Offers 1 switch visual access. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.




App	Device	About the App	Keystrokes	Switch Access
 <p>My Country \$9.99 (CA/US) \$10.49 (AU) £6.99 (UK)</p>	iPhone, iPad	A talking book about the United States of America. Key words are illustrated by SymbolStix. Users can also play bingo with symbols in the text, match symbols and play concentration with them. Part of the Click to Read series. Appropriate only to students in the USA.	Switches 1 – 4 can be any number key, ~ + 1, space or enter.	Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 <p>New Age Kurling \$2.99 (CA/AU/US) £1.99 (UK)</p>	iPad	A simple but addictive (for me!) version of Kurling for the iPad. Can be played with 1 or 2 players. Can be played at slow, medium or fast speed. Players need to aim their stones at the numbered target. The player with the most stones on the target at the end wins.	Switch – space or ~ + 1	Offers 1 switch visual access. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 <p>Rad Sounds \$4.99 (CA/US) \$5.49 (AU) £2.99 (UK)</p>	iPad	Rad Sounds is a cause/effect music program that is switch accessible. A customized playlist of music can be created by recording or singing new music to your iPad, importing from your iPad's music library or using some of the included songs. Switch access options include: Momentary mode, where the music plays for only the time you hold your finger on the iPad; Timed mode, where the music plays for a specified number of seconds upon each activation; or Advanced mode where the user chooses songs by cycling through their playlist.	Switch - ~ + 1	Cause and Effect. Can be accessed with RJ Cooper/Therapy Box and APPLICATOR switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.




App	Device	About the App	Keystrokes	Switch Access
 <p>Recycle It \$9.99 (CA/US) \$10.49 (AU) £6.99 (UK)</p>	iPhone, iPad	A talking book about recycling. Key words are illustrated by SymbolStix. Users can also play bingo with symbols in the text, match symbols and play concentration with them. Part of the Click to Read series. Appropriate only to students in the USA.	Switches 1 – 4 can be any number key, ~ + 1, space or enter.	Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 <p>Scan and Match \$19.99 (CA/US) \$20.99 (AU) £13.99 (UK)</p>	iPhone, iPad	A visual matching activity. The user is presented with a picture at the top of the screen, with a blank square next to it. At the bottom of the screen are two, three or four choices (selected in the settings), which can be all from the same category of vocabulary e.g. flowers or random across all the categories in the program. The user needs to select which one of the three choices is the same. Alternatively, the program can be set to a cause and effect mode where the correct picture goes up as soon as the switch or screen is selected.	Switches 1 – 4 can be any number key, ~ + 1, space or enter.	1, 2, or 3 switch visual scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 <p>Scan-A-Word \$4.99 (CA/US) \$5.49 (AU) £2.99 (UK) (Lite version available for free)</p>	iPad	Scan-A-Word is a single switch letter matching and spelling activity. A word (which may or may not be accompanied by a picture) appears on the screen. The first time it appears, a changing letter display appears below the first letter, and the student is encouraged to hit their switch when the matching letter appears. If the student selects the wrong letter, they start again. Once this task is completed, the word is repeated auditorily (or as a picture if one has been selected) and this time the student needs to select the correct letters to spell the word. Letters are consistently referred to by their letter names. Activity requires good auditory memory and isn't traditional scanning skills as the target constantly changes.	Switch - ~ + 1	Single switch visual scanning (not traditional scanning). Can be accessed with RJ Cooper/Therapy Box and APPLICATOR switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.




App	Device	About the App	Keystrokes	Switch Access
 <p>Scene&Heard \$49.99 (CA/US) \$51.99 (AU) £34.99 (UK) (Lite version available for free)</p>	iPhone, iPad	<p>Designed to construct communication books and pages with visual scene displays (VSDs) but can be used in a variety of creative ways as an educational or language tool. A range of pre-made displays are provided-or a user's own photos or Google images can be used at the VSD. Each VSD can be used as a standalone or can be linked to other scenes. Hybrid scenes can be created with a communication strip available on the bottom of the scene. Hotspots can play a recorded sound, pull up a symbol, launch a video or link to another scene. Once a scene has been made you are prompted to email it to yourself and others for back-up. Scenes can be browsed as images, laid out in a grid format, or scheduled to launch at specific times. No message window.</p>	<p>Switch 1 - ~ + 1 Switch 2 - ~ + 3</p>	<p>1 and 2 switch visual scanning. Can be accessed with RJ Cooper/Therapy Box and APPLICATOR switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.</p>
 <p>Sensory CineFx – Fun Photo Cinema Style with Hall of Mirrors, Psychedelic Colors and SpecialFX Free</p>	iPhone, iPad	<p>A photo effects app with 30 cinema style effects. A menu of four options can be accessed by switch 1, 2, 3 or 4. In each activity. Switch 1 flips the camera from forwards to backwards, switch 2 goes to the menu, switch 3 moves to a different effect and switch 4 takes a photo. If the user only has one switch site, the app could be used with a partner who could cycle through the graphics and the user could take photos of the ones they want (or vice versa). With 2 switches the user could have one switch (sending the number 3) which would move the image through the different effects (underwater, embedded in fire etc) and one switch (sending the number 4) would take a picture which is automatically added to the camera roll. If the user can access 3 or 4 switches then they can have access to the full features of each activity.</p>	<p>Switch 1 – 1 Switch 2 – 2 Switch 3 – 3 Switch 4 – 4</p>	<p>1 - 4 switch cause and effect. Can be accessed with the APPLICATOR switch interface or a Bluetooth keyboard. Also compatible with the Switch2Scan and SimplyWorks.</p>




App	Device	About the App	Keystrokes	Switch Access
 Sensory Electra Free	iPhone, iPad	<p>A simple cause and effect app showing an image of a light band across the screen. When the user presses any switch an audio track is played like an electrical shock is happening. The graphics remain brightly coloured on a black background. Switch 1 also changes the colours on the screen. Switch 2 moves the left side of the graphic up or down. Switch 3 moves the whole graphic up and down in a linear or non-linear fashion randomly. Switch 4 moves the right side of the graphic up or down. If the user only has one switch site, the app could be used with a partner who could cycle through the colours (switch 1) and the user could take move the image with switch 3 (or vice versa). With 2 switches the user could have one switch that would move the image through the different colours and one switch that would move the image around the screen. If the user can access 3 or 4 switches then they can have access to the more features of the app.</p>	Switch 1 – 1 Switch 2 – 2 Switch 3 – 3 Switch 4 – 4	1 - 4 switch cause and effect. Can be accessed with the APPLICATOR switch interface or a Bluetooth keyboard. Also compatible with the Switch2Scan and SimplyWorks.
 Sensory FotoFrez – A Fun Foto Booth with Hall of Mirrors, Psychedelic and Kaleidoscope Effects Free	iPhone, iPad	<p>A photo effects app with 30+ effects. A menu of four options can be accessed by switch 1, 2, 3 or 4. In each activity. Switch 1 creates a bubble ripple, switch 2 goes to the menu, switch 3 moves to a different effect and switch 4 takes a photo. If the user only has one switch site, the app could be used with a partner who could cycle through the graphics and the user could take photos of the ones they want (or vice versa). With 2 switches the user could have one switch (sending the number 3) which would move the image through the different effects (wavy or bubbly) and one switch (sending the number 4) would take a picture which is automatically added to the camera roll. If the user can access 3 or 4 switches then they can have access to the full features of each activity.</p>	Switch 1 – 1 Switch 2 – 2 Switch 3 – 3 Switch 4 – 4	1 - 4 switch cause and effect. Can be accessed with the APPLICATOR switch interface or a Bluetooth keyboard. Also compatible with the Switch2Scan and SimplyWorks.
 Sensory iMeba Free	iPhone, iPad	<p>Sensory iMeba starts with a brightly coloured blob on the screen. Each time the user presses the switch, the blob changes – it might elongate or double or triple. The graphics stay in the yellow and red colour range with a black background and some are accompanied by a sound effect. Each switch has a different range of effects. Could be used as a single or multi switch cause and effect app.</p>	Switch 1 – 1 Switch 2 – 2 Switch 3 – 3 Switch 4 – 4	1 - 4 switch cause and effect. Can be accessed with the APPLICATOR switch interface or a Bluetooth keyboard. Also compatible with the Switch2Scan and SimplyWorks.

App	Device	About the App	Keystrokes	Switch Access
 Sensory Just Touch Free	iPhone, iPad	A visual effects app with 30+ effects. A menu of four options can be accessed by switch 1, 2, 3 or 4 – option 1 uses the camera but options 2 – 4 are kaleidoscope style images. In each activity. Switch 1 creates changes the effect, switch 2 goes to the menu, switch 3 simulates a direct touch and switch 4 moves to the next effect. If the user only has one switch site, the app can easily be used to send keystroke 1, 2 or 4. If the user has more switch sites then more effects will be available.	Switch 1 – 1 Switch 2 – 2 Switch 3 – 3 Switch 4 – 4	1 - 4 switch cause and effect. Can be accessed with the APPLICATOR switch interface or a Bluetooth keyboard. Also compatible with the Switch2Scan and SimplyWorks.
 Sensory Magma Free	iPhone, iPad	Sensory iMeba starts with a brightly coloured blob on the screen. Each time the user presses the switch, the blob changes – it might elongate or double or triple. The graphics remain brightly coloured on a black background and are accompanied by a sound effect. Each switch has a different range of effects. Could be used as a single or multi switch cause and effect app.	Switch 1 – 1 Switch 2 – 2 Switch 3 – 3 Switch 4 – 4	1 - 4 switch cause and effect. Can be accessed with the APPLICATOR switch interface or a Bluetooth keyboard. Also compatible with the Switch2Scan and SimplyWorks.
 Sentence Match Chores: WHO is Doing WHAT \$3.99 (CA/US) \$4.49 (AU) £2.49 (UK)	iPad	A sentence and picture matching app. In the picture (or sentence) activity is taking place e.g. "He is clearing the table". In one mode a single picture can be matched to two, three or four sentences. For example, the picture "he is clearing the table" has the sentence options "He is clearing the table", "He is mowing the lawn" or "He is walking the dog" when three sentences are offered at the medium difficulty level. Alternatively a single sentence can be matched to two, three or four pictures. There are three difficulty levels. Speech feedback can be turned on or off. Sentences can be accompanied by Picture Communication Symbols (PCS) from Boardmaker or presented on their own. Scanning speed can be controlled in the settings.	In scanning mode - Switch 1 – Enter or Space or 1 or 2 or ~1 or ~2 In multi-switch mode Switch 1 – 1 Switch 2 – 2 Switch 3 – 3 Switch 4 – 4	1 switch scanning or 4 switches for selecting options 1 – 4. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.

App	Device	About the App	Keystrokes	Switch Access
 Shhhh! \$2.99 (CA/AU/US) £1.99 (UK)	iPad	A great app for someone learning to use two switches. Switch 1 opens the bedroom door to a noisy person e.g. a one-man band. Switch 2 makes the man in the bed sit up and get them to shhhh. No timing requirements.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Offers 2 switch visual access. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 Smarty Pants \$2.99 (CA/AU/US) £1.99 (UK)	iPad	This app from Inclusive Technology has large animated expressions and silly faces making this simple choice activity a fun way to explore feelings and encourage imitation. Ideal for children in first stages of learning to communicate.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Offers 1 and 2 switch visual access. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 Space Shooter \$2.99 (CA/AU/US) £1.99 (UK)	iPad	A simple single switch game with some timing requirements. Can be played at easy, medium or hard levels. Players need to press their switch when the target lines up with the spaceship. In the easy level, the time frame within which they need to press the switch goes for a couple of seconds, but if you increase the level of difficulty then this time frame decreases.	Switch – space or ~ + 1	Offers 1 switch visual access. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.

App	Device	About the App	Keystrokes	Switch Access
 <p>Special Stories \$13.99 (CA/US) \$14.99 (AU) £9.99 (UK)</p>	iPhone, iPad	An app for creating your own stories. User can create stories with photographs or pictures and then share them with others using Bluetooth, WiFi or email. Pages can be changed by touching the forward and back buttons on the screen or by pressing switch 1 to move forward. Switch 2 plays the sound for the page OR turns the page back, dependant on the settings you select.	Switch 1 – space Switch 2 – enter	Offers 1 and 2 switch visual access. Can be accessed with APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 <p>Special Words \$13.99 (CA/US) \$14.99 (AU) £9.99 (UK)</p>	iPhone, iPad	An app designed to have students matching words to pictures, words to words or pictures to pictures. The user needs to select a word (or picture) from a selection of 4. The setup of the app encourages a focus on nouns or words that are easily represented pictorially rather than high frequency words.	Switch 1 – space Switch 2 – enter	Offers 1 and 2 switch visual scanning. Can be accessed with APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 <p>Splat the clowns \$2.99 (CA/AU/US) £1.99 (UK)</p>	iPad	A simple single switch game with some timing requirements. Can be played at easy, medium or hard levels. Players need to press their switch within a given time frame to throw a pie or tomato (or other object) into the clown's face. In the easy level, the time frame within which they need to press the switch goes for a couple of seconds, but if you increase the level of difficulty then this time frame decreases.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Offers 1 switch visual access. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.



App	Device	About the App	Keystrokes	Switch Access
 <p>Switch Accessible Puzzles \$19.99 (CAN/US) \$20.99 (AU) £13.99 (UK)</p>	iPhone, iPad	A switch accessible puzzle app. User has the choice of number of puzzle pieces, with options from 2 to 5. Selection of pieces can be cause and effect, where the user simply presses the switch to have the piece go into the right spot – or scanning based where the user has to more actively complete the puzzle. App can also be used through direct access.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	1 switch cause and effect or 1 switch auto scan and 2 switch step scanning. Can be accessed with RJ Cooper/Therapy Box, Blue2 or APPLICATOR switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 <p>Switch Box Invaders \$2.99 (CA/AU/US) £1.49 (UK)</p>	iPad	A re-creation of the classic space invaders game. User has some control over speed, switch debounce time and there is a visual high contrast mode. A simple but engaging game.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 <p>SwitchKids \$19.99 (CAN/US) \$20.99 (AU) £13.99 (UK) (Lite version available for \$0.99)</p>	iPhone, iPad	Switch Kids contains three activities, Funny Sounds and Faces, Bubble Gum and Build-A-Kid. A scoreboard tracks the child's progress. Switch Kids teaches simple cause and effect, multiple-choice cause and effect, and a variation called "cause and effect with a goal".	Switches 1 – 4 can be any number key, ~ + 1, space or enter.	1, 2, 3 or 4 switch access with cause and effect activities. Can be accessed with RJ Cooper/Therapy Box, Blue2 or APPLICATOR switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.




App	Device	About the App	Keystrokes	Switch Access
 <p>SwitchTrainer \$16.99 (CA/US) \$17.99 (AU) £11.99 (UK) (Lite version available for free)</p>	iPad	SwitchTrainer is an app designed to provide the user to learn to use 1 or 2 switches, with a focus on moving towards 2 switches. Users can use external switches or the iPad screen divided into two halves. The app has 6 different activities: 1. a simple animation activity where switch 1 triggers the animals (e.g. a monkey playing bongos) and switch 2 moves between animations; 2. An activity using the same animations, but where the user has a choice of 2 and can select which one to display; 3. A playground where the user can move a monkey around and then have them do different activities; 4. An activity where the user chooses the parts of a costume to dress a mouse; 5. An activity where the user can conduct a small band, choosing which instrument plays; 5. Hide and seek where the user has to find one of the characters hiding behind some objects. The app offers control over timing, number of switches and switch debounce.	Switch 1 – space Switch 2 – enter	Offers 1 and 2 switch visual scanning. Can be accessed with APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. iPad screen can also be used as 1 or 2 switches. Also compatible with Switch2Scan and SimplyWorks.
 <p>Train Tracker \$2.99 (CA/AU/US) £1.99 (UK)</p>	iPad	Train tracker requires the user to select a train – and then build a track for it to travel around. 7 different track pieces can be selected and added to build a track with curves and an intersection. Once the track is built, the train can travel around it and the user can take a snapshot of their track.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Offers 1 and 2 switch visual scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 <p>Treasure Island \$2.99 (CA/AU/US) £1.99 (UK)</p>	iPad	A map appears on the screen with 9 segments. The user can scan or use switch 1 to move around the map, trying to guess where the treasure is stored. The map is presented differently each time and the treasure is stored in a different location.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.




Specialbites now has some iPad switch accessible videos and slideshows - <http://www.specialbites.com/index.html> and see the link on the left hand side.









Switch Accessible Apps (AAC)




App	Device	About the App	Keystrokes	Switch Access
 <p>Alexicom AAC Free Subscription to Alexicom online USD\$40 per month or \$408 per annum</p>	iPhone, iPad	Vocabulary is organized in a list view. Vocabulary can be organized categories, alphabetically, in recency or "in popular" order. User can add an unlimited number of vocabulary items using their own photos and recording their own speech. No message window for sentence construction. A number of page sets available from Alexicom Online. Text-to-speech includes word prediction.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.
 <p>Alexicom Elements Adult Home (Female) \$19.99 (CAN/US) \$20.99 (AU) £13.99 (UK)</p>	iPad	A multi-level dynamic display 16 cell page set with a mixture of situation specific or activity based vocabulary and core vocabulary. A message window allows the user to put together utterances using single words, phrases and typing. Use of the Apple keyboard is supported by additional word prediction. This version has vocabulary/topics/photographs aimed at a adult female with an emphasis on vocabulary for use at home. Easy editing of both existing items and new items with option of adding new pages as needed. App interfaces with Twitter and includes data tracking and ability to backup to iCloud.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.




App	Device	About the App	Keystrokes	Switch Access
 Alexicom Elements Adult Home (Male) \$19.99 (CAN/US) \$20.99 (AU) £13.99 (UK)	iPad	<p>A multi-level dynamic display 16 cell page set with a mixture of situation specific or activity based vocabulary and core vocabulary. A message window allows the user to put together utterances using single words, phrases and typing. Use of the Apple keyboard is supported by additional word prediction. This version has vocabulary/topics/photographs aimed at a adult male with an emphasis on vocabulary for use at home. Easy editing of both existing items and new items with option of adding new pages as needed. App interfaces with Twitter and includes data tracking and ability to backup to iCloud.</p>	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICator or Blue2 switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.
 Alexicom Elements Child Home (F) SymbolStix \$19.99 (CAN/US) \$20.99 (AU) £13.99 (UK)	iPad	<p>A multi-level dynamic display 16 cell page set with a mixture of situation specific or activity based vocabulary and core vocabulary. A message window allows the user to put together utterances using single words, phrases and typing. Use of the Apple keyboard is supported by additional word prediction. This version has vocabulary/topics/photographs aimed at a girl and includes SymbolStix symbols. Easy editing of both existing items and new items with option of adding new pages as needed. App interfaces with Twitter and includes data tracking and ability to backup to iCloud.</p>	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICator or Blue2 switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.
 Alexicom Elements Child Home (Female) \$19.99 (CAN/US) \$20.99 (AU) £13.99 (UK)	iPad	<p>A multi-level dynamic display 16 cell page set with a mixture of situation specific or activity based vocabulary and core vocabulary. A message window allows the user to put together utterances using single words, phrases and typing. Use of the Apple keyboard is supported by additional word prediction. This version has vocabulary/topics/photographs aimed at a girl. Easy editing of both existing items and new items with option of adding new pages as needed. App interfaces with Twitter and includes data tracking and ability to backup to iCloud.</p>	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICator or Blue2 switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.



App	Device	About the App	Keystrokes	Switch Access
 <p>Alexicom Elements Child Home (M) SymbolStix \$19.99 (CAN/US) \$20.99 (AU) £13.99 (UK)</p>	iPad	A multi-level dynamic display 16 cell page set with a mixture of situation specific or activity based vocabulary and core vocabulary. A message window allows the user to put together utterances using single words, phrases and typing. Use of the Apple keyboard is supported by additional word prediction. This version has vocabulary/topics/photographs aimed at a boy and includes SymbolStix symbols. Easy editing of both existing items and new items with option of adding new pages as needed. App interfaces with Twitter and includes data tracking and ability to backup to iCloud.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.
 <p>Alexicom Elements Child Home (Male) \$19.99 (CAN/US) \$20.99 (AU) £13.99 (UK)</p>	iPad	A multi-level dynamic display 16 cell page set with a mixture of situation specific or activity based vocabulary and core vocabulary. A message window allows the user to put together utterances using single words, phrases and typing. Use of the Apple keyboard is supported by additional word prediction. This version has vocabulary/topics/photographs aimed at a boy. Easy editing of both existing items and new items with option of adding new pages as needed. App interfaces with Twitter and includes data tracking and ability to backup to iCloud.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.
 <p>Alexicom Elements Child Pre (F) SymbolStix \$19.99 (CAN/US) \$20.99 (AU) £13.99 (UK)</p>		A multi-level dynamic display 9 cell page set with a mixture of situation specific or activity based vocabulary and core vocabulary. A message window allows the user to put together utterances using single words, phrases and typing. Use of the Apple keyboard is supported by additional word prediction. This version has vocabulary/topics/photographs aimed at a preschool girl and includes SymbolStix symbols. Easy editing of both existing items and new items with option of adding new pages as needed. App interfaces with Twitter and includes data tracking and ability to backup to iCloud.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.




App	Device	About the App	Keystrokes	Switch Access
 <p>Alexicom Elements Child Pre (Female) \$19.99 (CAN/US) \$20.99 (AU) £13.99 (UK)</p>	iPad	<p>A multi-level dynamic display 9 cell page set with a mixture of situation specific or activity based vocabulary and core vocabulary. A message window allows the user to put together utterances using single words, phrases and typing. Use of the Apple keyboard is supported by additional word prediction. This version has vocabulary/topics/photographs aimed at a preschool girl. Easy editing of both existing items and new items with option of adding new pages as needed. App interfaces with Twitter and includes data tracking and ability to backup to iCloud.</p>	<p>Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3</p>	<p>Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.</p>
 <p>Alexicom Elements Child Pre (M) SymbolStix \$19.99 (CAN/US) \$20.99 (AU) £13.99 (UK)</p>	iPad	<p>A multi-level dynamic display 9 cell page set with a mixture of situation specific or activity based vocabulary and core vocabulary. A message window allows the user to put together utterances using single words, phrases and typing. Use of the Apple keyboard is supported by additional word prediction. This version has vocabulary/topics/photographs aimed at a preschool boy and includes SymbolStix symbols. Easy editing of both existing items and new items with option of adding new pages as needed. App interfaces with Twitter and includes data tracking and ability to backup to iCloud.</p>	<p>Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3</p>	<p>Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.</p>
 <p>Alexicom Elements Child Pre (Male) \$19.99 (CAN/US) \$20.99 (AU) £13.99 (UK)</p>	iPad	<p>A multi-level dynamic display 9 cell page set with a mixture of situation specific or activity based vocabulary and core vocabulary. A message window allows the user to put together utterances using single words, phrases and typing. Use of the Apple keyboard is supported by additional word prediction. This version has vocabulary/topics/photographs aimed at a preschool boy. Easy editing of both existing items and new items with option of adding new pages as needed. App interfaces with Twitter and includes data tracking and ability to backup to iCloud.</p>	<p>Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3</p>	<p>Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.</p>




App	Device	About the App	Keystrokes	Switch Access
 <p>Alexicom Elements Core 280 (Female) \$9.99 (CAN/US) \$10.49 (AU) £6.99 (UK)</p>	iPad	Core vocabulary based app. Main page has 54 items, mostly high frequency words. These are generally accompanied by a symbol which is the word itself on a coloured background. Would be helpful to have more symbols on the pages as it can be very slow to locate words. Words are coloured using Fitzgerald key e.g. yellow for nouns etc. Easy to edit and customise. Message window allows construction of utterances using core vocabulary and typing. App interfaces with Twitter and includes data tracking and ability to backup to iCloud.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.
 <p>Alexicom Elements Core 280 (Male) \$9.99 (CAN/US) \$10.49 (AU) £6.99 (UK)</p>	iPad	Core vocabulary based app. Main page has 54 items, mostly high frequency words. These are generally accompanied by a symbol which is the word itself on a coloured background. Would be helpful to have more symbols on the pages as it can be very slow to locate words. Words are coloured using Fitzgerald key e.g. yellow for nouns etc. Easy to edit and customise. Message window allows construction of utterances using core vocabulary and typing. App interfaces with Twitter and includes data tracking and ability to backup to iCloud.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.
 <p>Alexicom Elements Story Maker \$9.99 (CAN/US) \$10.49 (AU) £6.99 (UK)</p>	iPad	A collection of story templates designed to enable the easy sharing of stories. Each story is 8 cells or 8 pages, with 425 story pages in total available plus the option for adding more pages as needed. Each item can have text, photograph/image and speech. Very easy to use and customise – making construction of a story while out and about a definite option. User can also add additional comments using text to speech and the Apple onscreen keyboard at any time and the app interfaces with Twitter. Typing is supported by word prediction and has inbuilt data tracking. Stories can be backed up to iCloud.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.



App	Device	About the App	Keystrokes	Switch Access
 <p>Alexicom Elements Teen Home (Female) \$19.99 (CAN/US) \$20.99 (AU) £13.99 (UK)</p>	iPad	<p>A multi-level dynamic display 16 cell page set with a mixture of situation specific or activity based vocabulary and core vocabulary. A message window allows the user to put together utterances using single words, phrases and typing. Use of the Apple keyboard is supported by additional word prediction. This version has vocabulary/topics/photographs aimed at a teen female. Easy editing of both existing items and new items with option of adding new pages as needed. App interfaces with Twitter and includes data tracking and ability to backup to iCloud.</p>	<p>Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3</p>	<p>Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.</p>
 <p>Alexicom Elements Teen Home (Male) \$19.99 (CAN/US) \$20.99 (AU) £13.99 (UK)</p>	iPad	<p>A multi-level dynamic display 16 cell page set with a mixture of situation specific or activity based vocabulary and core vocabulary. A message window allows the user to put together utterances using single words, phrases and typing. Use of the Apple keyboard is supported by additional word prediction. This version has vocabulary/topics/photographs aimed at a teen male with an emphasis on vocabulary for use at home. Easy editing of both existing items and new items with option of adding new pages as needed. App interfaces with Twitter and includes data tracking and ability to backup to iCloud.</p>	<p>Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3</p>	<p>Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.</p>
 <p>Alexicom Elements UK Child Home (Female) \$19.99 (CAN/US) \$20.99 (AU) £13.99 (UK)</p>	iPad	<p>A multi-level dynamic display 16 cell page set with a mixture of situation specific or activity based vocabulary and core vocabulary. A message window allows the user to put together utterances using single words, phrases and typing. Use of the Apple keyboard is supported by additional word prediction. This version has vocabulary/topics/photographs aimed at a girl. Easy editing of both existing items and new items with option of adding new pages as needed. This version includes UK spelling and accent. App interfaces with Twitter and includes data tracking and ability to backup to iCloud.</p>	<p>Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3</p>	<p>Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.</p>



App	Device	About the App	Keystrokes	Switch Access
 <p>Alexicom Elements UK Child Home (Male) \$19.99 (CAN/US) \$20.99 (AU) £13.99 (UK)</p>	iPad	A multi-level dynamic display 16 cell page set with a mixture of situation specific or activity based vocabulary and core vocabulary. A message window allows the user to put together utterances using single words, phrases and typing. Use of the Apple keyboard is supported by additional word prediction. This version has vocabulary/topics/photographs aimed at a boy. Easy editing of both existing items and new items with option of adding new pages as needed. This version includes UK spelling and accent. App interfaces with Twitter and includes data tracking and ability to backup to iCloud.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.
 <p>Alexicom Elements UK Teen Home (Female) \$19.99 (CAN/US) \$20.99 (AU) £13.99 (UK)</p>	iPad	A multi-level dynamic display 16 cell page set with a mixture of situation specific or activity based vocabulary and core vocabulary. A message window allows the user to put together utterances using single words, phrases and typing. Use of the Apple keyboard is supported by additional word prediction. This version has vocabulary/topics/photographs aimed at a teen female. Easy editing of both existing items and new items with option of adding new pages as needed. This version includes UK spelling and accent. App interfaces with Twitter and includes data tracking and ability to backup to iCloud.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.
 <p>Alexicom Elements UK Teen Home (Male) \$19.99 (CAN/US) \$20.99 (AU) £13.99 (UK)</p>	iPad	A multi-level dynamic display 16 cell page set with a mixture of situation specific or activity based vocabulary and core vocabulary. A message window allows the user to put together utterances using single words, phrases and typing. Use of the Apple keyboard is supported by additional word prediction. This version has vocabulary/topics/photographs aimed at a teen male. Easy editing of both existing items and new items with option of adding new pages as needed. This version includes UK spelling and accent. App interfaces with Twitter and includes data tracking and ability to backup to iCloud.	Switch 1 – space or ~ + 1 Switch 2 – enter or ~ + 3	Visual and auditory 1 and 2 switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.



App	Device	About the App	Keystrokes	Switch Access
 <p>GoTalk Now for iPad \$79.99 (CAN/US) \$84.99 (AU) £54.99 (UK) (Free version available)</p>	iPad	<p>AAC App which allows the user to create a number of pages with 1, 4, 9, 16 or 25 items per page, or a visual scene with multiple buttons. Single pages or a multipage book with page links can be created. Pictures or videos can be displayed on each cell. A core vocabulary row of 4 items can be accessed easily from every page. Visual scene displays can be created using hot spots. Create multiple users/books, share the books via wifi, create a pdf of pages and email them, use the GoTalk Now template library or email links to templates you have created. Backup to a computer or to DropBox. Easy to use and create unlimited pages. Has message window for sentence construction on express pages - but this isn't available on all page types. A number of pre-made pages can be downloaded - these have been made by a range of different people.</p>	<p>Switch 1 - ~ + 1 Switch 2 - ~ + 3 (But need to have settings set to RJ Cooper switch interface)</p>	<p>1 or 2 switch visual or auditory scanning. Can be accessed with RJ Cooper/Therapy Box or APPLICATOR switch interfaces or a Bluetooth keyboard. The Attainment Switch App on an iPhone or second iPad. Auditory scanning includes option for private auditory scan cues using one audio channel. Also compatible with Switch2Scan and SimplyWorks.</p>
 <p>MyTalk Tools Mobile \$99.99 (CAN/US) \$109.99 (AU) £69.99 (UK) MyTalk Workspace US\$75 per annum (Lite version available for free)</p>	iPhone, iPad	<p>Comes with limited pages. MyTalk Workspace needs to be purchased to allow you to design your own pages, including access to an Image Library, Web Authoring and collaboration tools. No message window for sentence construction. A cell can be set up to pull-up the apple keyboard for text-to-speech, but this cannot be incorporated with symbol images to form a sentence. Video modeling using the camera in the iPad2.</p>	<p>Switch - ~ + 1 or space</p>	<p>Visual and auditory single switch scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.</p>




App	Device	About the App	Keystrokes	Switch Access
 PhotoVOCA Free Free	iPhone, iPad	PhotoVOCA AAC is a multi-level communication app. It comes with a couple of example vocabulary sets which are category based. Easy to use and program and allows flexible creation of a range of pages. No message window.	Switches 1 – 4 can be any number key, ~ + 1, space or enter.	Can be accessed with RJ Cooper/Therapy Box, Blue2 or APPLICATOR switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch with 1 or 2 finger taps to emulate switches. Also compatible with Switch2Scan and SimplyWorks.
 PicBoard \$7.99 (CAN/US) \$8.49 (AU) £5.49 (UK)	iPhone, iPad	Can design individual communication boards with anything from 1 large cell to 20 cells. Easy to use and includes a nice "how to" video. This is the non-speaking version of TalkBoard. Boards can be printed and easily shared with others. Visual timers can be created on communication pages. Can import Picture Communication Symbols from Boardmaker using the free PicCutter tool for PC and Mac.	Switch 1 - ~ + 1 Switch 2 - ~ + 3	1 and 2 switch, visual and auditory scanning. Can be accessed with RJ Cooper/Therapy Box or APPLICATOR switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 Proloquo2Go \$219.99 (CAN/US) \$229.99 (AU) £149.99 (UK)	iPhone, iPad	Full featured AAC app with two default language options. A comprehensive core vocabulary arrangement gives options for using a grid from 9 to 64 buttons on each page. A "basic communication" option is aimed at one to two word users and can be arranged in between 9 and 36 cells per page. Every cell can be customized and the user has access to SymbolStix picture symbols as well as the option of importing their own photos. A message window allows the user to construct sentences by stringing words together. A keyboard can be used to type for text-to-speech and typing can inserted into the message window and added to a sentence. Word prediction is available to support typing. User editing allows the easy creation of extra pages, which can be linked into the existing page structure or used to create a new user from scratch. Multiple users are supported.	1 or 2 switches can be set to a keystroke of the user's choice from a range of keystrokes	Can be accessed with RJ Cooper/Therapy Box, Blue2 or APPLICATOR switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch with either the full screen functioning as 1 switch or 2 halves functioning as 2 switches. Also compatible with Switch2Scan and SimplyWorks.




App	Device	About the App	Keystrokes	Switch Access
 SayHi! AAC Free	iPhone, iPad	Simple AAC App. User can create multiple pages which each link off a central menu. 9 cells can be displayed on each page. Each cell can display text, a picture and have recorded speech. Cells can also be easily setup to play an audio file installed on the iPad. The App can be installed on three devices on the same Wi-Fi network and the second two devices can be used for alternative access, functioning as if they are the two switches in two switch scanning.	NA	2 switch visual scanning using two other iPhones or iPod touch devices running on the same Wi-Fi network. Also compatible with Switch2Scan and SimplyWorks.
 Scan To Pictures \$4.99 (CAN/US) \$5.49 (AU) £2.99 (UK) Lite version available for free.	iPad	Designed as an AAC Training App, Scan to Pictures offers a series of options on a page - between 2 and 9. Once the user selects an item, the app tells them which item they chose and can then play an optional consequence or reinforcer. The same 2 - 9 items are presented until edited. If the user doesn't make a selection, the App plays images and sounds to tempt them to interact with it.	Switch - ~ + 1	Single switch visual scanning and auditory scanning. Can be accessed with RJ Cooper/Therapy Box and APPLICATOR switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.
 Scene&Heard \$49.99 (CAN/US) \$51.99 (AU) £34.99 (UK)	iPhone, iPad	Designed to construct communication books and pages with visual scene displays (VSDs) but can be used in a variety of creative ways as an educational or language tool. A range of pre-made displays are provided-or a user's own photos or Google images can be used at the VSD. Each VSD can be used as a standalone or can be linked to other scenes. Hybrid scenes can be created with a communication strip available on the bottom of the scene. Hotspots can play a recorded sound, pull up a symbol, launch a video or link to another scene. Once a scene has been made you are prompted to email it to yourself and others for back-up. Scenes can be browsed as images, laid out in a grid format, or scheduled to launch at specific times. No message window.	Switch 1 - ~ + 1 Switch 2 - ~ + 3	1 and 2 switch visual scanning. Can be accessed with RJ Cooper/Therapy Box or APPLICATOR switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.




App	Device	About the App	Keystrokes	Switch Access
 <p>Sounding Board Free</p>	iPhone	Sounding Board is an App from AbleNet that allows you to create communication displays from 1 to 9 message locations using AbleNet symbols or your own photos; can link displays to other displays; allows auditory scanning in addition to single and dual switch scanning; record messages of virtually any length for a symbol/photo.	Switch 1 – space, 1 or ~ + 1 Switch 2 – enter, 3 or ~ + 3	1 and 2 switch visual and auditory scanning. Can be accessed with RJ Cooper/Therapy Box, Blue2 or APPLICator switch interfaces or a Bluetooth keyboard. Also compatible with Switch2Scan and SimplyWorks.
 <p>TalkBoard \$16.99 (CAN/US) \$17.99 (AU) £11.99 (UK)</p>	iPhone, iPad	Can design individual communication boards with anything from 1 large cell to 20 cells. Easy to use and includes a nice "how to" video. This is the speaking version of PicBoard. Boards can be printed and easily shared with others. Visual timers can be created on communication pages. Can import Picture Communication Symbols from Boardmaker using the free PicCutter tool for PC and Mac.	Switch 1 - ~ + 1 Switch 2 - ~ + 3	1 and 2 switch, visual and auditory scanning. Can be accessed with RJ Cooper/Therapy Box or APPLICator switch interfaces or a Bluetooth keyboard. Auditory scanning includes option for private auditory scan cues using one audio channel. Also compatible with Switch2Scan and SimplyWorks.


App	Device	About the App	Keystrokes	Switch Access
 <p>TalkRocket Go \$99.99 (CAN/US) \$109.99 (AU) £69.99 (UK)</p>	iPhone, iPad	<p>User can customise the pre-existing vocabulary and add their own at www.myvoiceaac.com or on the device it is being used on. Any changes made to the user's page set online automatically sync with device, as long as you are logged in as same user in both places. Comes pre-installed with conversation starter based "books" with phrases like "I want". User may need to select down through multiple levels before app speaks e.g. "I want" then "to go" then "to the bathroom" and then finally speak "to floss", although individual items can be spoken as you go by swiping over them. Vocabulary sets on different topics also available online to be added e.g. basketball. Can use GPS feature in iPhone/iPad to mark a location and add a photo to it. Once a place has been added on the device, the user can create specific vocabulary for that location. User can toggle between text-to-speech mode and symbol/picture mode but cannot combine items from both to construct sentences.</p>	NA	<p>1 switch visual scanning – automatic or step scanning. Can be accessed with RJ Cooper/Therapy Box or APPLICATOR switch interfaces or a Bluetooth keyboard. Use iPad screen as a switch. Also compatible with Switch2Scan and SimplyWorks.</p>
 <p>TalkTablet NEO \$79.99 (CAN/US) \$84.99 (AU) £54.99 (UK)</p>	iPhone, iPad	<p>Very simple to use app for creating a multi-level communication system. Comes with a pre-programmed page sets for 16, 24 and 32 cells aimed at emerging AAC users. Existing page sets can also be customised or a new page set can be created from scratch with unlimited pages and between 1 – 256 cells per page. Multiple users can be easily created with different settings and pages. A message window allows the user to string together words to speak a sentence or phrase. A custom keyboard can be created to enable the user to input text into the same message window. The keyboard mode on the toolbar uses the system keyboard, and provides good text-to-speech, but this operates with a different message window. Includes iVerbalChat for communication wirelessly with other TalkTablet users over any distance.</p>	<p>Switch 1 – space, 1 or ~ + 1 Switch 2 – enter, 3 or ~ + 3</p>	<p>1 and 2 switch linear or row/column visual and auditory scanning. Can be accessed with RJ Cooper/Therapy Box or APPLICATOR switch interfaces or a Bluetooth keyboard. Can use iPad screen as a switch. Also compatible with Switch2Scan and SimplyWorks.</p>

App	Device	About the App	Keystrokes	Switch Access
 TalkTablet UK \$89.99 (CAN/US) \$94.99 (AU) £59.99 (UK)	iPhone, iPad	Very simple to use app for creating a multi-level communication system. Comes with a pre-programmed page sets for 16, 24 and 32 cells aimed at emerging AAC users. Existing page sets can also be customised or a new page set can be created from scratch with unlimited pages and between 1 – 256 cells per page. Multiple users can be easily created with different settings and pages. A message window allows the user to string together words to speak a sentence or phrase. A custom keyboard can be created to enable the user to input text into the same message window. The keyboard mode on the toolbar uses the system keyboard, and provides good text-to-speech, but this operates with a different message window. Includes iVerbalChat for communication wirelessly with other TalkTablet users over any distance.	Switch 1 – space, 1 or ~ + 1 Switch 2 – enter, 3 or ~ + 3	1 and 2 switch linear or row/column visual and auditory scanning. Can be accessed with RJ Cooper/Therapy Box or APPLICATOR switch interfaces or a Bluetooth keyboard. Can use iPad screen as a switch. Also compatible with Switch2Scan and SimplyWorks.
 TalkTablet US \$89.99 (CAN/US) \$94.99 (AU) £59.99 (UK)	iPhone, iPad	Very simple to use app for creating a multi-level communication system. Comes with a pre-programmed page sets for 16, 24 and 32 cells aimed at emerging AAC users. Existing page sets can also be customised or a new page set can be created from scratch with unlimited pages and between 1 – 256 cells per page. Multiple users can be easily created with different settings and pages. A message window allows the user to string together words to speak a sentence or phrase. A custom keyboard can be created to enable the user to input text into the same message window. The keyboard mode on the toolbar uses the system keyboard, and provides good text-to-speech, but this operates with a different message window. Includes iVerbalChat for communication wirelessly with other TalkTablet users over any distance.	Switch 1 – space, 1 or ~ + 1 Switch 2 – enter, 3 or ~ + 3	1 and 2 switch linear or row/column visual and auditory scanning. Can be accessed with RJ Cooper/Therapy Box or APPLICATOR switch interfaces or a Bluetooth keyboard. Can use iPad screen as a switch. Also compatible with Switch2Scan and SimplyWorks.

App	Device	About the App	Keystrokes	Switch Access
 <p>TapSpeak Button Standard for iPad \$29.99 (CAN/US) \$31.99 (AU) £20.99 (UK)</p>	iPad	Turns your iPad into a single message communication device/switch. Allows the user to record a number of messages and store them – which can then be pulled up again latter. Has an option of tap-feedback when the button is pressed.	Switch - ~ + 1 or space	Cause and effect. Can be accessed with RJ Cooper/Therapy Box, Blue2 or APPLICator switch interfaces or a Bluetooth keyboard. Can use iPad screen as a switch. Also compatible with Switch2Scan and SimplyWorks.
 <p>TapSpeak Button Plus for iPad \$49.99 (CAN/US) \$51.99 (AU) £34.99 (UK)</p>	iPad	Turns your iPad into a single message communication device/switch. Allows the user to record a number of messages and store them- which can then be pulled up later. Has an option of tap-feedback when the button is pressed. TapSpeak Button Plus includes the Picture Communication Symbol library from Mayer Johnson.	Switch - ~ + 1 or space	Cause and effect. Can be accessed with RJ Cooper/Therapy Box, Blue 2 or APPLICator switch interfaces or a Bluetooth keyboard. Can use iPad screen as a switch. Also compatible with Switch2Scan and SimplyWorks.
 <p>TapSpeak Choice for iPad \$149.99 (CAN/US) \$159.99 (US) £99.99 (UK)</p>	iPad	AAC app for building multi-level or single board communication systems. Users need to create a library with vocabulary in it - and then select items from the library to be added to a page. Can have between 1 and 56 cells per board and can arrange vocabulary with gaps to keep items in specific locations. Keyboard and text only buttons can be combined with symbol buttons to produce sentences in the message window.	Switch 1 - ~ + 1 or space Switch 2 - ~ + 3 or enter	1 or 2 switch visual and auditory scanning. Can be accessed with RJ Cooper/Therapy Box, APPLICator or Blue2 switch interfaces or a Bluetooth keyboard. Can use iPad screen as a switch. Also compatible with Switch2Scan and SimplyWorks.




App	Device	About the App	Keystrokes	Switch Access
 <p>TapSpeak Sequence Standard \$29.99 (CAN/US) \$31.99 (AU) £20.99 (UK)</p>	iPhone, iPad	Turns your iPad into a sequenced single message communication device/switch. Allows the user to record a number of sequences and store them – which can then be used again later. Can change how long the “switch” needs to be pressed before it activates. Includes ability to backup and restore the sequence database.	Switch - ~ + 1	Cause and effect. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Can use iPad screen as a switch. Also compatible with Switch2Scan and SimplyWorks.
 <p>TapSpeak Sequence Plus \$49.99 (CAN/US) \$51.99 (AU) £34.99 (UK)</p>	iPhone, iPad	Turns your iPad into a sequenced single message communication device/switch. Allows the user to record a number of sequences and store them – which can then be used again later. Can change how long the “switch” needs to be pressed before it activates. Includes ability to backup and restore the sequence database. TapSpeak Sequence Plus includes the Picture Communication Symbol library from Mayer Johnson.	Switch - ~ + 1	Cause and effect. Can be accessed with RJ Cooper/Therapy Box, APPLICATOR or Blue2 switch interfaces or a Bluetooth keyboard. Can use iPad screen as a switch. Also compatible with Switch2Scan and SimplyWorks.
 <p>TapToTalk Free Tap to Talk Designer \$99.95USD for one year or \$179.95 outright</p>	iPhone, iPad	Comes with a starter album with categories “go, play, hungry, drink, emotional, bathroom”. Can make your own albums on-line with the TapToTalk Designer. No message window for sentence construction. TapToTalk Uploader lets the user add photos, pictures direct from the iPad/iPhone to their TapToTalk Designer Account.	Switch - ~ + 1	Single switch visual scanning. Can be accessed with RJ Cooper/Therapy Box or APPLICATOR switch interfaces or a Bluetooth keyboard. Can use iPad screen as a switch. Also compatible with Switch2Scan and SimplyWorks.


App	Device	About the App	Keystrokes	Switch Access
 <p>Predictable \$159.99 (CAN, US) \$169.99 (AU) £109.99 (UK)</p>	iPhone, iPad	Full featured text-to-speech app, with word prediction with the ability to type text on an app-specific QWERTY keyboard and email it as well as speaking it. On the iPhone version the text can be sent as a txt message. You can also use the text you have typed to do a status updated in Facebook or to Tweet. It allows you to save what you have typed as a phrase under a range of categories, and also pull up a history of previous conversations and has an option of hand writing input.	Switch 1 - ~ + 1 or space Switch 2 - ~ + 3 or enter	1 or 2 switch visual scanning with some auditory scanning (but not complete). Can be accessed with RJ Cooper/Therapy Box or APPLICATOR switch interfaces or a Bluetooth keyboard. Can use iPad screen as a switch. Also compatible with Switch2Scan and SimplyWorks.
 <p>Voice4U \$59.99 (CAN, US) \$64.99 (AU) £39.99 (UK)</p>	iPhone iPad	Vocabulary is organised in a list view. Vocabulary can be organised in categories, alphabetically, in recency or in "popular" order. User can add an unlimited number of vocabulary items using their own photos and recording their own speech. No message window for sentence construction.	Switch - ~ + 1	Single switch visual scanning. Can be accessed with RJ Cooper/Therapy Box and APPLICATOR switch interfaces or a Bluetooth keyboard. iPad screen can be used as a switch. Also compatible with Switch2Scan and SimplyWorks.
 <p>iMyVoice - AAC \$39.99 (CAN, US) \$41.99 (AU) £27.99 (UK)</p>	iPad	Simple multi-level AAC app that lets you create a range of pages that link to each other. Each page can have between 1 and 100 cells on it. Button can contain both images and text. A button can speak, link to another page or open a web page, movies or play music. No message window for creating sentences. Easy to use. Small pre-made page set which is mostly based around needs and wants.	iPad screen only	Single switch visual scanning using iPad screen as switch

App	Device	About the App	Keystrokes	Switch Access
 <p>UcanTALK \$2.99 (CA, US, AU) £1.99 (UK)</p>	iPad	A text based AAC app. User can create a number of categories and store colour coded items under each category. Items are labelled with text only and have recordings attached to them. A text section allows the user to pull up the system keyboard and type a phrase for their communication partner to read (no voice output in this section).	-	Not compatible with existing switch interfaces. Company seems to be developing their own.



Switch Interfaces Switch Accessible Apps for iPad/iPhone

Picture	Switch Interface and Links	Cost
 <p>(Has switch port for second switch to be used in two switch scanning)</p>	RJ Cooper Bluetooth Super Switch	\$149USD from RJ Cooper
 <p>(Can use inbuilt switches, or has two switch jacks for external switches)</p>	Blue 2 Bluetooth Switch	\$149USD from AbleNet \$219AUD from Spectronics
 <p>(Has two switch ports for external switches)</p>	AirTurn BT105 Also sold as: SwitchBox from TherapyBox RJ Cooper Bluetooth Switch Interface	\$99USD from RJ Cooper \$155AUD from Zyteq or \$189 from Spectronics

 <p>(Has four ports for external switches. Can be used to control music as well as Apps.)</p>	<p>APPLICATOR</p>	<p>\$165USD from Inclusive</p> <p>\$280AUD from Zytec or \$329AUD from Spectronics</p>
--	-----------------------------------	--



NB – As a cheap, one switch option, you can also store a single keystroke (e.g. space) in [USB Button](#) and then connect it to the iPad with a USB adaptor to perform the functions of a single switch. (Thanks to John Gale for this tip).

If you are interested in using a PC switch interface to control these apps, check out this [informative video from Bill Zeigler](#).



Switch Interfaces for VoiceOver and VoiceOver compatible apps

These are a group of switch interfaces which can be used with VoiceOver to provide switch access to the iOS. Some of these also provide switch access to the switch accessible apps. A Bluetooth keyboard can also be used to control VoiceOver.

Picture	Switch Interface and Links	Cost
 <p>(Has two switch ports and an extra adaptor can be used to provide ports for up to 6 switches)</p>	<p>Tecla Shield is available from Komodo. Jane has reviewed this at http://www.janefarrall.com/blog/2012/05/29/tecla-shield-for-iphone-and-ipad-adding-greater-switch-access-to-ios/</p>	<p>\$289CAD from Komodo</p> <p>\$349.99AUD from Spectronics</p>
 <p>(Has four switch ports)</p>	<p>Switch2Scan is available from Pretorian and Inclusive. Jane has reviewed this at http://www.janefarrall.com/blog/2012/06/29/switch2scan-switch-access-to-apps-ios-and-more/ . The Switch2Scan can also be used to access all of the switch accessible apps listed above by changing to a different mode within the Switch2Scan.</p>	<p>\$180USD from Inclusive</p> <p>\$275 from Zytec or \$379AUD from Spectronics</p>

 <p>(No switch ports – works with the wireless products in the SimplyWorks range)</p>	<p>SimplyWorks for iPad is part of the wireless Simply Works range of access products from Pretorian. This can be used to access switch accessible apps plus gives VoiceOver control to iOD. Jane has reviewed this at http://www.janefarrall.com/blog/2012/09/30/simplyworks-for-ipad-part-1-switch-access/. The real strength of this interface is with switch accessible apps and media or access to iOS with the joystick http://www.janefarrall.com/blog/2012/10/04/simplyworks-for-ipad-part-2-joystick-and-keyboard-access/.</p>	<p>\$195USD from Inclusive</p> <p>\$300 from Zytec or \$409AUD from Spectronics</p>
	<p>Zybox for iOS VoiceOver interface. Available with a 30 pin port or Lightning connector. This interface has been reviewed by Jane at http://www.janefarrall.com/blog/2013/06/08/zybox-for-ios-updated-and-upgraded/</p>	<p>\$345USD from Zygo</p> <p>Also available from Zytec</p>
	<p>iPad VoiceOver Controller from RJ Cooper – a multiple button controller option.</p>	<p>\$199 from RJ Cooper</p>

Upcoming Interfaces

[Perrero](#) by RSL Steeper

[J-Pad](#) from Pretorian

[Connect](#) from AbleNet